

# Kompatibilitätsliste

## Elektroporationsküvetten

Elektroporator	1 mm Küvette Kurze Elektrode 75-01-SE	2 mm Küvette Kurze Elektrode 75-02-SE	1 mm Küvette Lange Elektrode 75-01-LE	2 mm Küvette Lange Elektrode 75-02-LE	4 mm Küvette Lange Elektrode 75-04-LE
<b>Amaxa</b>					
Nucleofector	-	-	-	-	X
<b>Bio-Rad</b>					
Genepulser	X	X	X	X	X
Genepulser 2	X	X	X	X	X
Genepulser Xcell	-	-	X	X	X
Micropulser	X	X	X	X	X
E.coli Pulsar	X	X	X	X	X
<b>BTX</b>					
ECM 2001	X	X	X	X	X
ECM 399	X	X	X	X	X
ECM 630	X	X	X	X	X
ECM 830	X	X	X	X	X
Gemini SC2	X	X	X	X	X
Gemini X2	X	X	X	X	X
Agile Pulse	X	X	X	X	X
<b>Eppendorf</b>					
Electroporator 2510	X	X	X	X	X
Multiporator 4308	X	X	X	X	X
Eporator 4309	X	X	X	X	X